




Two Striper®

Snowpost® & Snowlight®

Two Striper® *Snowpost*® & *Snowlight*® posts are elastic modulus with flexural strength close to dentin.

Made of zircon-rich glass fiber embedded in a resin matrix.

- Superior fatigue resistance
- Radiopaque and light transmitting
- Autoclavable up to 132°C / 269.6°F
- Bio-compatible and corrosion free
- Anatomical shape and 5* different diameters (0.8mm* – 1.6mm) *Two Striper® *Snowlight*® posts only



Mechanical Properties		
	Two Striper® Snowpost®	Two Striper® Snowlight®
Elastic Modulus	45 GPa	52 GPa
Shear Strength	38 MPa	64 MPa

Two Striper® Snowpost® Kits

C810001	Basic Kit : 50 assorted posts (1.0, 1.2, 1.4, 1.6) and 4 matching finishing drills
C810002	Intro Kit: 20 assorted posts (1.0, 1.2, 1.4, 1.6) and 4 matching finishing drills

Two Striper® Snowpost® Refills

10 Posts/Pack		10 posts/pack with matching finishing drill	
C811001	Ø 1.0	C811005	Ø 1.0
C811002	Ø 1.2	C811006	Ø 1.2
C811003	Ø 1.4	C811007	Ø 1.4
C811004	Ø 1.6	C811008	Ø 1.6

Two Striper® Snowlight® Kits

C830001	Basic Kit : 50 assorted posts (1.0, 1.2, 1.4, 1.6) and 4 matching finishing drills
C830002	Intro Kit: 20 assorted posts (1.0, 1.2, 1.4, 1.6) and 4 matching finishing drills

Two Striper® Snowlight® Refills

10 Posts/Pack		10 posts/pack with matching finishing drill	
C831009	Ø 0.8		
C831002	Ø 1.2	C831006	Ø 1.0
C831003	Ø 1.4	C831007	Ø 1.2
C811004	Ø 1.6	C831008	Ø 1.4
C811004	Ø 1.6	C831009	Ø 1.6